


T3 PROSTHESIS

Cement-retained Solutions

T3

Angled Abutments



Angled abutments offer various solutions for different prosthetic platform diameters and gingiva heights in the cases of cemented angled restorations.

3.8 platform diameter has the options of 10°/20°/30° angles. 5.0 and 6.0 platform diameters have 20° angle option.

The diameter appropriate for planning can be chosen from the schema.

Abutments are produced of Gr5 titanium material.


Abutments are packaged with abutment screw.



A	3.8	3.8	3.8
B	0	0	0
Code	T3 4510-10	T3 4510-20	T3 4510-30



A	5.0	5.0
B	1.0	2.0
Code	T3 4511-20	T3 4512-20



A	6.0	6.0	
B	1.0	2.0	Abutment Screw
Code	T3 5011-20	T3 5012-20	T0 2001

