


T3 PROSTHESIS

Overdenture Solutions

T3

Ball Abutment



Ball abutment solutions are available for overdenture prosthetics with different gingiva heights.

Ball abutments are packed come with a silicone ring and matrix.

Abutments are produced of Gr5 titanium material.



A	5.0	5.0
B	1.5	3.0
Code	T3 5030	T3 5031

Instruments

Impression Post T0 32615	Analog T0 4521

Closed Metal Housing T0 5207	Mid-Hard Plastic T0 4500-16	Soft Plastic T0 4500-7	Metal Ring&Silicon T0 5200	Silikon Ring T0 4500-3

Metal Cup&Plastic T0 5201	Metal Cup Plastic T0 4500-4	Impression Post Silicon T0 4500-18

