


T3 PROSTHESIS

Auxiliary Prosthesis Parts

T3

Gingiva Formers




NucleOSS gingiva formers are designed in different height and sizes to achieve the design soft tissue shape.

Soft tissue shape enable to take impression easy, fast and comfortably, for both open and closed tray techniques.

It's made of Titanium Gr5 material. It is placed and taken out by hex driver.

A	5.0	5.0	6.0	6.0
B	1.0	4.0	1.0	4.0
Kod	T3 4516	T3 4517	T3 5016	T3 5017

Temporary Gingiva Formers



Suitable for all diameters.

