


T3 PROSTHESIS

Custom Solutions

T3

Universal Ti-Base Casting Abutments



Universal Ti-base casting abutments are offered for dentist and laboratories who design custom abutments with casting technique.

Abutments are produced of titanium Gr5 material.

Hex for single tooth and No-hex for multiple teeth restorations are present.

It is offered as a white burn-out cap and two screws(T3 16001).

Multiple Teeth Restorations

NO-HEX



A

5.0

B

0.9

Code

T3 5080

Single Teeth Restorations

HEX



5.0

0.9

T3 5081

Instruments



Casting Cap

T3 32907



Casting Cap

T0 32907



Abutment Screw


T3 16001

Impression Posts



Impression Post(Closed)

T3 4940-1



Open Impression(HEX)

T3 4940



Open Impression (NO-HEX)

T3 4941

Analogs



3.4 3.8 4.2

T3 4520



5.0

T3 5020

