


Gingiva Formers



The use of platform switch gingiva formers needs caution. They are designed to create a perfect relationship between the soft tissue and abutment. Therefore selecting the appropriate length and diameter of the gingiva former is essential. This will not only create a good relationship with the gum but avoid your patients to experience pain during try in sessions.

A recommended matching gingiva former & abutment table is below.

All diameters of gingiva formers can be used on all implant diameters if the case has special requirements.

Platform Switch

A	4.3	4.3	4.3	5.0	5.0	5.0
B	2.0	4.0	6.0	2.0	4.0	6.0
Code	T6 NR501	T6 NR502	T6 NR503	T6 SD504	T6 SD505	T6 SD506

A	6.0	6.0	6.0
B	2.0	4.0	6.0
Code	T6 WD507	T6 WD508	T6 WD509


Cylindrical Form

A	4.3	4.3	4.3	5.0	5.0	5.0
B	2.0	4.0	6.0	2.0	4.0	6.0
Code	T6 NR512	T6 NR513	T6 NR514	T6 SD515	T6 SD516	T6 SD517

A	6.0	6.0	6.0
B	2.0	4.0	6.0
Code	T6 WD518	T6 WD519	T6 WD520

Temporary Abutment

6.0
T6 32034



*The material is PEEK and it is suitable for shaping.