


Impression Posts



Impression posts are used for easy, fast and comfortable impression making.

There are two different designs for closed & open tray techniques.

Hex & No-Hex designs are available.

Over-Implant Impression Posts

	HEX	NO-HEX	HEX	NO-HEX
A	4.3	4.3	4.3	4.3
B	2.0	2.0	2.0	2.0
Code	T6 32611	T6 32612	T6 32604	T6 32605

Impression Cap
T0 32912

Over-Abutment Impression Posts

NO-HEX	NO-HEX	HEX	HEX
Impression Post (Close)	Impression Post (Open)	Impression Post (Close)	Impression Post (Open)
T0 32607	T0 32608	T0 32609	T0 32610

Connection Screws					
	T6 16000	T6 16002	T6 16006	T6 16003	T6 16004

Analogs	NO-HEX	HEX
	T6 32201	T0 32202
		T0 32203

