


T6 PROsthESIS

Cement-retained Solutions



Standard Straight Abutments



Provides varying gingival heights & diameters for cemented restorations.

Information about the abutments can be seen in below table.

Abutment packages include impression cap, abutment screw to be used for the patient and abutment screw to be used in the laboratory (yellow).

They are made of Gr5 Titanium material.

A	4.3	4.3	4.3	4.3	4.3
B	1.0	2.0	3.0	4.0	5.0
Code	T6 NR040	T6 NR041	T6 NR042	T6 NR043	T6 NR044

A	5.0	5.0	5.0	5.0	5.0
B	1.0	2.0	3.0	4.0	5.0
Code	T6 SD045	T6 SD046	T6 SD047	T6 SD048	T6 SD049

A	6.0	6.0	6.0
B	1.5	3.0	4.5
Code	T6 WD050	T6 WD051	T6 WD052

Impression Cap			
	4.3	5.0	6.0
	T0 32909	T0 32910	T0 32911

Abutment Screw	Laboratory Screw
T6 16000	T6 16007

