








Gingiva Formers




NucleOSS gingiva formers are designed as various height and sizes for gingiva to be in the best shape during healing process.


Produced of Titanium Gr5 material. It is placed and is taken out by hex driver.

			
A	4.2	4.2	4.2
B	2.0	4.0	6.0
Code	T44052	T44054	T44056

			
A	5.0	5.0	5.0
B	2.0	4.0	6.0
Code	T44852	T44854	T44856

			
A	6.0	6.0	6.0
B	2.0	4.0	6.0
Code	T46052	T46054	T46056

Temporary Gingiva Former


6.0
T432034

