


# Tpure PROSTHESIS

## Screw-retained Solutions



### Multi Unit Abutment | No-Hex



Screw-retained abutments offer accurate and permanent solutions for different prosthetic platform diameters and gingiva heights in the cases of screw-retained restorations.

0°, 17° ve 30° angled options are present.

Abutments are produced of titanium Gr5 material.

**Straight** Multi unit abutments are offered with a **white** burn out cap and two occlusal screws.

**Angled** Multi Unit abutments are offered with a **white** burn out cap, a multi unit connection screw and two occlusal screws.



A	5.0	5.0	5.0	5.0	5.0	5.0	5.0
B/C	1.5	3.0	5.0	1.5/3.1	2.4/4.0	1.25/4.0	2.0/4.75
Code	T45045	T45046	T45047	T45541-17	T45542-17	T45543-30	T45544-30

Instruments (No-Hex)

Comfort Cap	Scan Post	Ti-Base Abutment	Titanium Cap	Burn out Cap
T0 32906	T0 32033	T0 32811	T0 32913	T0 32901

Impression Post (Open)	Multi Unit Analog	Implant Connection Screw	Abutment Occlusal Screw	(Open) Impression Screw
T0 32608	T0 32202	T42004	T0 16003	T0 16004

Impression Post (Closed)	Impression Cap
T0 32607	T0 32912