


### Straight Abutments



Straight abutments offer various solutions for different prosthetic platform diameters and gingiva heights in the cases of cemented restorations.

The diameter appropriate for planning can be chosen from the schema.


Abutments are produced of titanium Gr5 material.

Abutments are packaged with abutment screw.

A	4.2	4.2	4.2	4.2
B	1.0	2.0	3.0	4.0
Code	T44001	T44002	T44003	T44005

A	5.0	5.0	5.0	5.0
B	1.0	2.0	3.0	4.0
Code	T44801	T44802	T44803	T44805

A	6.0	6.0	6.0	6.0
B	1.0	2.0	3.0	4.0
Code	T46001	T46002	T46003	T46005



Temporary Straight Abutments

5.0	5.0
HEX	NO-HEX
T432022	T432023

Impression Caps

4.3	5.0	6.0		
T0 32909	T0 32910	T0 32911	Impression Cap	Abutment Screw
			T0 32912	T0 2001