


Universal Ti-Base Casting Abutments

Universal Ti-base casting abutments are offered for dentist and laboratories who design custom abutments with casting technique.

Abutments are produced of titanium Gr5 material.

Hex for single tooth and No-hex for multiple teeth restorations are present.

It is offered as a white burn-out cap and two screws(T4 16001).



Multiple Teeth Restorations

	NO-HEX	NO-HEX
A	4.0	4.5
B	0.9	0.9
Code	T44080	T44580

Single Teeth Restorations

	HEX	HEX
A	4.0	4.5
B	0.9	0.9
Code	T44081	T44581

Instruments

Casting Cap	Casting Cap	Abutment Screw
T432907	T0 32907	T416001

Impression Posts

5.0	5.0	5.0
Impression Post(Closed)	Open Impression(HEX)	Open Impression (NO-HEX)
T45540-1	T45540	T45541

Analogs

3.4 3.8 4.2	5.0
T44020	T44820

